Disappearing Ink

Purpose

To illustrate the effect of pH change on indicator color.

Materials

- Thymolphthalein
- Ethanol
- Sodium hydroxide solution, NaOH (0.1 M)
- 100-mL beaker or bottle

Safety

- Read the SDS sheets for all chemicals before using them.
- Wear safety glasses and gloves.
- Sodium hydroxide is caustic; handle with care.

Procedure

- 1. Dissolve a small amount of thymolphthalein in ~10 mL of ethanol.
- 2. Slowly add 40 mL of water while stirring.
- 3. Add aqueous sodium hydroxide solution, 0.1 M, dropwise, until solution turns dark blue.
- 4. The resulting solution can be used as disappearing ink.

Results

• The ink disappears over time.

Follow-up Teaching Notes

- Thymolphthalein is blue above pH 10 and colorless below pH 9.
- The ink solution becomes colorless as it absorbs carbon dioxide from the air

$$2 \text{ NaOH} + \text{CO}_2 \rightarrow \text{Na}_2\text{CO}_3 + \text{H}_2\text{O}$$

Connections

• acid/base indicators, pH.

Extension

• Two common invisible inks are:

AgNO_{3 (aq)}, which can be developed by light CuSO_{4 (aq)}, which can be developed by ammonia

Disposal/Clean-up

- The paper can be disposed of in the garbage.
- Extra ink can be placed in a sealed and properly labeled container for reuse.

Ward's Science Tel: (866) 260-0501

105